WKt 6451

Wine storage cabinet **GrandCru**



LIEBHERR

Freestanding

WKt 6451

Wine storage cabinet **GrandCru**



← 75 →



Energy consump. in 365 days/24 h	143 / 0.391 kWh
Max. No. of Bordeaux bottles 0.75 Ltr	3121
Frequency / voltage	50 Hz / 220-240 V~
Connection Rating	1.0 A
Door	Full door, terra
Side walls	terra
Exterior dimensions in cm (h / w / d)	193 / 74.7 / 75.9
Controls	
Type of Control	MagicEye on front, key electronic
Temperature display	Digital temperature display
Alarm in case of malfunction	visual and acoustic
Door alarm	acoustic
Child proof function	yes
Features	
Temperature zones	1
Temperature range	+5 °C to +20 °C
Adjustable cooling circuits	1
Forced air system	Convection cooling, fresh air supply via FreshAir activated
	charcoal filter
Interior light	LED interior lighting
Handle	Aluminium handle with integrated opening mechanism
Design	SwingLine-Design
Shelf material	Adjustable beech wood shelves
Storage shelves	7
Door hinges	Reversible door hinging
Lock	fitted
EAN-Nr.	9005382238677

LIEBHERR

 $^{^1\}text{These}$ models are only bound for storing wine. Maximum DIN norm bordeaux bottles capacity (0.75 l) as per norm NF H 35-124 (H: 300.5 mm, Ø 76.1 mm)

Freestanding WKt 6451

Features





Door alarm

The audible door alarm alerts the user if the door has been open for longer than 60 seconds for safe protection of the food.



Child lock

The child lock is programmed so as to prevent the appliance being inadvertently switched off. A symbol in the MagicEye indicates when the child lock is switched on.



Reversible door hinges

The appliances are factory fitted with door hinges on the right. The reversible door hinges permit optimum use of the appliance wherever it is located.



Wine storage cabinets

Wine storage cabinets provide similar conditions to wine cellars. A constant temperature prevails throughout the interior in these appliances. It can be set to between +5 °C and +20 °C as required. Wine storage cabinets are ideal for the long-term storage and maturing of wines. Depending on the set temperature, these appliances can however also store more ample stocks of wine at serving temperature.



Activated charcoal filter

Environmental odours can affect wines stored for longer periods. In the wine cabinets, perfect air quality is ensured through an easily-replaceable FreshAir activated charcoal filter.





Low-vibration storage

Even slight vibrations can disturb the wine maturing process and prevent the tannins from settling. Specially developed, exceptionally low-vibration compressors in all the wine cabinets mean that the wine is optimally stored.

LIEBHERR

Technical drawing





Products and their details are staged for advertising purposes. Whilst every effort is made to ensure all specifications and information in this flyer are accurate, Liebherr reserves the right to make changes without notification. Where actual appliance measurements, specifications, and ratings are crucial, we recommend that the actual appliance installation instructions that accompany each appliance be checked. Liebherr will not be liable for any damage. (01.12.2022)





