


# WKb 4112

Wine storage cabinet

Barrique



# LIEBHERR

Freestanding

# WKb 4112

Wine storage cabinet  
Barrique

breite\_66.pdf



Wine

Energy consump. in 365 days/24 h	178 kWh
Max. No. of Bordeaux bottles 0.75 Ltr	168 <sup>1</sup>
Frequency / voltage	60 Hz / 115 V~
Connection Rating	1.5 A
Door	Glass door, frame black
Side walls	black
Exterior dimensions in cm (h / w / d)	164.5 / 65.7 / 66.3
<b>Controls</b>	
Type of Control	MagicEye on front, key electronic
	Digital temperature display
Temperature display	Digital temperature display
Alarm in case of malfunction	visual and acoustic
Door alarm	acoustic
<b>Features</b>	
Temperature zones	1
Temperature range	+5 °C to +20 °C
Adjustable cooling circuits	1
Forced air system	Convection cooling, fresh air supply via FreshAir activated charcoal filter
Interior light	LED interior lighting
Handle	Ergonomically designed slimline handle, black
Design	HardLine-Design
Shelf material	Adjustable beech wood shelves
Storage shelves	6
Pull-out shelves on telescopic rails	1
Door hinges	Reversible door hinging
Lock	fitted
EAN-Nr.	9550000011380

<sup>1</sup> These models are only bound for storing wine. Maximum DIN norm bordeaux bottles capacity (0.75 l) as per norm NF H 35-124 (H: 300.5 mm, Ø 76.1 mm)

## Features



### Door alarm

The audible door alarm alerts the user if the door has been open for longer than 60 seconds for safe protection of the food.



### Reversible door hinges

It is possible to change the door hinge on the appliance. The reversible door hinges permit optimum use of the appliance wherever it is located.



### Wine storage cabinets

Wine storage cabinets provide similar conditions to wine cellars. A constant temperature prevails throughout the interior in these appliances. It can be set to between +5 °C and +20 °C as required. Wine storage cabinets are ideal for the long-term storage and maturing of wines. Depending on the set temperature, these appliances can however also store more ample stocks of wine at serving temperature. In the GrandCru and Vinothek series, Liebherr offers multi-temperature wine cabinets




### Activated charcoal filter

Environmental odours can affect wines stored for longer periods. In the wine cabinets, perfect air quality is ensured through an easily-replaceable FreshAir activated charcoal filter.



### Presentation shelves

The presentation shelf sets off fine wines to perfection and keeps bottles that have already been opened at the optimum temperature. Up to 6 wine bottles can be placed upright on the flexible shelf. Wine bottles can be stored lying down either behind or at the side.




### Low-vibration storage

Even slight vibrations can disturb the wine maturing process and prevent the tannins from settling. Specially developed, exceptionally low-vibration compressors in all the wine cabinets mean that the wine is optimally stored.

## Technical drawing

---



Products and their details are staged for advertising purposes. Whilst every effort is made to ensure all specifications and information in this flyer are accurate, Liebherr reserves the right to make changes without notification. Where actual appliance measurements, specifications, and ratings are crucial, we recommend that the actual appliance installation instructions that accompany each appliance be checked. Liebherr will not be liable for any damage. (03.03.2023)

