



Quality Wines Deserve
Quality Storage

Quality, Design and Innovation

Wkt 6451

LIEBHERR



Electronic LCD Display: Intuitive navigation of the button control system enables tailored storage options.



Wooden shelves: Liebherr's natural untreated beechwood shelves are ideal for storing Bordeaux bottles. Store bottles in opposite directions to make full use of the capacity.



Solid door: A solid door wine cabinet mimics the conditions of an underground wine cellar, making it perfect for maturing wines by creating a dark, undisturbed environment with zero exposure to UV light and heat fluctuations.

The **WKt 6451** wine storage offers similar conditions to a wine cellar. Temperature can be set to exactly what you require, within the range of +5°C to +20°C, and is maintained constantly and uniformly throughout the interior.

Liebherr wine cabinets are perfect for storing wines for long periods of time to allow them to mature to perfection. Depending on what temperature the appliances are set to, they can also be used for keeping white wine or red wine at exactly the right temperature for serving. Our integrated activated charcoal filter¹ absorbs odours from the ambient air and protects the taste of your fine wines; best changed once every 12 months.


WKt 6451

GrandCru

Wine

Sound rating dB(A) *	41
Total capacity litre	625
Max No. of Bordeaux bottles (0.75l) **	312
Number of shelves	7
Shelf material	Beech wood
Electronic control	MagicEye on front, buttons
Temperature range °C	+5°C to +20°C
Temperature zones	1
Temperature alarm	Visual and audible
Interior light	LED light bulb
Door alarm	Audible
Safety thermostat	Yes
Activated charcoal filter	Yes
Door type	Full door, terra
Integrated lock	Yes
SoftSystem	No
Door hinges	Right, reversible
Product dimensions (H/W/D) cm	193 / 74.7 / 75.9
Gross / Net weight kg	100 / 87
Voltage V/Hz	220-240 / 50

Installation Diagram



Activated Charcoal Filter
Part no.: 988129100

* Noise output - according to EN 60704-3

** These models are only bound for storing wine. Max. DIN norm bordeaux bottle capacity (0.75 l) as per norm NF H 35-124 (H:300.5 mm, Ø 76.1mm)

¹ Available accessory: Activated Charcoal Filter (Part no. 988129100),

Note: Whilst every effort is made to ensure all specifications and information in this flyer are accurate, Liebherr reserves the right to make changes without notification. Where actual appliance measurements, specifications, and ratings are crucial, we recommend that the actual appliance installation instructions that accompany each appliance be checked. Liebherr will not be liable for any loss.