

Savour the Difference
with a Liebherr

Quality, Design and Innovation



UWTgb 1682

LIEBHERR



TipOpen: When lightly tapped, TipOpen technology gently opens the glass door by 7cm to facilitate further opening. If accidentally brushed open, the door automatically closes after 3 seconds (if it has not been pulled open further).



Wooden shelf: The practical wooden pull-out shelves on telescopic rails enable smooth shelf sliding and are ideal for storing Bordeaux bottles. Store your wine bottles in opposite directions to make full use of the capacity.



Electronic LCD Display: Intuitive navigation of touch electronic control system enables tailored storage options. The digital temperature display can be viewed via the UV tinted glass door.


The dual-zone **UWTgb 1682** wine cabinet, with its unique black glass frame, is designed to be a highlight in your home. Our favourite feature is the innovative TipOpen mechanism that allows you the convenience of opening the appliance door with just one gentle tap! Our natural, handcrafted, untreated beechwood shelves release no chemical aromas and reduce vibrations so that your wine can develop in peace.

This elegant undercounter model will fit perfectly into your kitchen island or bar counter. With two independent temperature zones, you will be able to store different varieties of wines simultaneously, each at its optimum drinking or serving temperature. An activated charcoal filter¹ protects stored wines against odour; best changed once every 12 months.

UWTgb 1682
Vinidor

LED	Wine	Tip Open	Soft System	Telescopic Rails
-----	------	----------	-------------	------------------

Sound rating dB(A) *	36
Total net capacity litre	94
Max No. of Bordeaux bottles (0.75l) **	34
Number of shelves	5
Shelves on telescopic rails	3
Shelf material	Beech wood shelves
Electronic control	LCD display, touch electronic
Temperature range °C	+5°C to +20°C
Temperature zones	2
Temperature alarm	Visual and audible
Interior light	LED
Door alarm	Audible
Safety thermostat	Yes
Activated charcoal filter	Yes
Door type	Flushed Black Glass Door
SoftSystem	Yes
Door hinges	Right, reversible
Product dimensions (H/W/D) cm	81.8 / 59.7 / 57.7
Cut out dimensions (H/W/D) cm	82 - 87 / 60 / min. 58
Gross / Net weight kg	52 / 48
Voltage V/Hz	220-240 / 50



Activated Charcoal Filter
Part no.: 988128900

* Noise output - according to EN 60704-3

** These models are only bound for storing wine. Max. DIN norm bordeaux bottle capacity (0.75 l) as per norm NF H 35-124 (H:300.5 mm, Ø 76.1mm)

¹ Available accessory: Activated Charcoal Filter (Part no. 988128900)

Note: Whilst every effort is made to ensure all specifications and information in this flyer are accurate, Liebherr reserves the right to make changes without notification. Where actual appliance measurements, specifications, and ratings are crucial, we recommend that the actual appliance installation instructions that accompany each appliance be checked. Liebherr will not be liable for any loss.