


The Perfect Environment for Your Wine

Quality, Design and Innovation



Wkr 4211

LIEBHERR



Optimum humidity: Liebherr wine cabinets guarantee an optimum humidity and this ensured the corks are supple in place whilst still allowing the wine to breathe in its maturing process.



Wooden shelves: Liebherr's natural untreated beechwood shelves are ideal for storing Bordeaux bottles. Stored bottles in opposite directions to make full use of the capacity.



Solid door: A solid door wine cabinet mimics the conditions of an underground wine cellar, making it perfect for maturing wines by creating a dark, undisturbed environment with zero exposure to UV light and heat fluctuations.


The **WKr 4211** solid door wine cabinet mimics the conditions of an underground wine cellar perfect for maturing wines by creating a dark, undisturbed environment with zero exposure to UV light and heat fluctuations.

All Liebherr wine cabinets ensure a constant interior temperature that can be adjusted between +5°C and +20°C. Our handcrafted, untreated beechwood shelves release no chemical aromas and reduce vibrations so that your wine can develop in peace.

An activated charcoal filter¹ protects stored wines against odour and maintains air quality in the wine cabinet; best changed every 12 months.

WKr 4211	Wine
Vinothek	

Sound rating dB(A) *	39
Total net capacity litre	377
Max No. of Bordeaux bottles (0.75l) **	200
Number of shelves	6
Shelf material	Beech wood
Electronic control	LCD display, buttons
Temperature range °C	+5°C to +20°C
Temperature zones	1
Temperature alarm	Visual and audible
Interior light	Light bulb
Door alarm	Audible
Safety thermostat	Yes
Activated charcoal filter	Yes
Door type	Full door, Bordeaux Red
Door hinges	Right, reversible
Product dimensions (H/W/D) cm	165 / 60 / 73.9
Gross / Net weight kg	74 / 67
Voltage V/Hz	220-240 / 50



Activated Charcoal Filter
Part no.: 988129100

* Noise output - according to EN 60704-3

** These models are only bound for storing wine. Max. DIN norm bordeaux bottle capacity (0.75 l) as per norm NF H 35-124 (H:300.5 mm, Ø 76.1mm)

¹ Available accessory: Activated Charcoal Filter (Part no. 988129100)

Note: Whilst every effort is made to ensure all specifications and information in this flyer are accurate, Liebherr reserves the right to make changes without notification. Where actual appliance measurements, specifications, and ratings are crucial, we recommend that the actual appliance installation instructions that accompany each appliance be checked. Liebherr will not be liable for any loss.