



ENERGY STAR CERTIFIED

Lab Grade Refrigerators and Freezers

Liebherr - SRFfg 3501 : SRFfg 3501

Specifications

ENERGY STAR Unique ID:	2668046
ENERGY STAR Partner:	Liebherr-Hausgerate Ochsenhausen GmbH
Brand Name:	Liebherr
Model Name:	SRFfg 3501
Model Number:	SRFfg 3501
Product Type:	High Performance Refrigerator
Total Volume (cubic feet):	11.19
Lab Grade Refrigerator or Freezer Energy Consumption (kWh/day):	0.63
Peak Temperature Variance (°C):	4.78
Product Form Factor:	Upright
Total Number of Doors:	1
Number of glass doors:	0
Number of solid doors:	1
Door Opening Orientation:	Swinging
Number of Drawers or Shelves:	4
Number of Internal Ports and/or Access Holes:	1
Depth (in):	26.14
Height:	66.3
Width:	23.5
Defrost Type:	Manual
Refrigerant with GWP:	R-600a (GWP:3)
Installed Accessories:	n/a
Date Available on Market:	1706763600000
Date Certified:	1695182400000
Markets:	United States, Canada
ENERGY STAR Certified:	Yes

Captured On:
06/10/2024