




# Wine Storage Sophistication For Your Home

Quality, Design and Innovation

# LIEBHERR



**Optimum humidity:** Liebherr wine cabinets guarantee an optimum humidity and this ensured the corks are supple in place whilst still allowing the wine to breathe in its maturing process.



**Integrated lock:** Nicely concealed into the side of the door, this sturdy lock protects your fine wine collection.



**Wooden shelves:** The wooden shelves are ideal for storing Bordeaux bottles. They can be stored in opposite directions to make full use of the capacity.

The **WKb 4611** wine cabinet mimics the conditions of an underground wine cellar. Its solid door creates a dark, undisturbed environment with zero exposure to UV light and heat fluctuations – perfect for maturing wines. This unit ensures a constant interior temperature, which can be adjusted between +5 °C and +20 °C. Wines are always stored safely on our untreated wooden shelves while the activated charcoal filter provides optimal air quality. The ergonomic slimline door handle with integrated opening mechanism makes the door easy to open with very little effort.

**WKb 4611**


Barrique

LED

Wine

Energy consumption (kWh/y)	175
Sound rating dB(A)	40
Total capacity litre	336
Max No. of Bordeaux bottles (0.75l) *	195
Number of shelves	7
Shelves on telescopic rails	0
Shelf material	Beech wood shelves
Presentation shelf	0
Electronic control	LCD display, key electronic
Temperature range °C	+5°C to +20°C
Temperature zones	1
Temperature alarm	Visual and audible
Interior light	Light bulb
Door alarm	Audible
Safety thermostat	Yes
Activated charcoal filter	Yes
Door type	Full door, black
Integrated lock	Yes
SoftSystem	No
Door hinges	Right, reversible
Product dimensions (H/W/D) cm	184.2 / 65.7 / 66.3
Gross / Net weight kg	78.5 / 68.5
Voltage V/Hz	220-240 / 50

**Installation Diagram**



Activated Charcoal Filter  
Part no.: 988129100

\* Noise output - according to EN 60704-3

\*\* These models are only bound for storing wine. Max. DIN norm bordeaux bottle capacity (0.75 l) as per norm NF H 35-124 (H:300.5 mm, Ø 76.1 mm)