




Perfect Display For  
All Your Favourites

Quality, Design and Innovation

**LIEBHERR**



**Optimum humidity:** Liebherr wine cabinets guarantee an optimum humidity and this ensured the corks are supple in place whilst still allowing the wine to breathe in its maturing process.




**LED lighting:** The tinted glass door provides not only a perfect UV protection, but also an optimum glance at wine. The integrated LED lighting ensures a pleasant illumination of the complete interior.



**Presentation shelf:** The presentation shelf, equipped with wooden angles, is ideal for presenting wine and storing open bottles safely. Up to 6 bottles can be set upright on the flexible shelf. In addition further bottles can be stored lying down either behind or at the side.

The interior of Liebherr's **WKb 4612** wine cabinet from the Barrique series has a constant temperature that can be set as needed between +5 °C to +20 °C. Wines are stored safely at all times on untreated wooden shelves while the activated charcoal filter provides optimal air quality. With the capacity to hold 195 Bordeaux bottles, this spacious model is ideal for the long-term storage and maturation of wines. And thanks to its coated glass door, your prized collection is always well protected from UV radiation. Efficient LED lighting ensures that the interior is illuminated evenly and pleasantly. Add-on accessories<sup>1</sup> are also available for added flexibility.

WKb 4612	
Barrique	
LED	Wine
Energy consumption (kWh/y)	204
Sound rating dB(A) *	40
Total capacity litre	336
Max No. of Bordeaux bottles (0.75l) **	195
Number of shelves	7
Shelves on telescopic rails	0
Shelf material	Beech wood shelves
Presentation shelf	1
Electronic control	LCD display, key electronic
Temperature range °C	+5°C to +20°C
Temperature zones	1
Temperature alarm	Visual and audible
Interior light	LED
Door alarm	Audible
Safety thermostat	Yes
Activated charcoal filter	Yes
Door type	Tinted glazed door, black frame
Integrated lock	Yes
SoftSystem	No
Door hinges	Right, reversible
Product dimensions (H/W/D) cm	184.2 / 65.7 / 66.3
Gross / Net weight kg	96 / 86
Voltage V/Hz	220-240 / 50



Activated Charcoal Filter  
Part no.: 988129100

\* Noise output - according to EN 60704-3

\*\* These models are only bound for storing wine. Max. DIN norm bordeaux bottle capacity (0.75 l) as per norm NF H 35-124 (H:300.5 mm, Ø 76.1mm)

<sup>1</sup> Optional accessories: Shelf to store wine upright (Part no. 936674700), Retrofitting Telescopic Shelf (Part no. 945575000), Activated Charcoal Filter (Part no: 988129100)