

Created For True
Wine Connoisseurs

Quality, Design and Innovation



LIEBHERR



Electronic LCD display: Intuitive navigation of touch electronic control system enables tailored storage options. The digital temperature display can be viewed via the UV tinted glass door.



Presentation shelf: The practical shelf of our Vinidor range fulfils both bottle storage and presentation functions, while telescopic rails enable smooth shelf sliding. Bottles can be stored horizontally at the rear.



LED lighting: The light in the cooling compartments can be dimmed to give a pleasing ambient light inside the appliance.


The Liebherr **WTes 5872** multi-temperature wine cabinet satisfies your thirst for the smartest, most elegant ways to store your precious wines by providing ideal conditions for them to mature successfully. This is achieved with consistent air quality, temperature, and humidity control alongside a quiet, vibration-free storage. Additionally an activated charcoal filter protects your finest collection against odour; best changed once every 12 months.

With three adjustable temperature zones and controlled by an innovative LCD display, you can simultaneously store red wine, white wine, and champagne, all at their optimum serving temperature. Impress your guests with what is a true wine storage unit!

WTes 5872
Vinidor

- LED
- Wine
- Soft System
- Telescopic Rails

Energy consumption (kWh/y)	245
Sound rating dB(A) *	38
Total capacity litre	511
Max No. of Bordeaux bottles (0.75l) **	178
Number of shelves	13
Shelves on telescopic rails	10
Shelf material	Beech wood shelves
Presentation shelf	1
Electronic control	LCD display, touch electronic
Temperature range °C	+5°C to +20°C
Temperature zones	3
Temperature alarm	Visual and audible
Interior light	LED (Dimmable)
Door alarm	Audible
Safety thermostat	Yes
Activated charcoal filter	Yes
Door type	Glazed door with stainless steel frame
Integrated lock	Yes
SoftSystem	Yes
Door hinges	Right, reversible
Product dimensions (H/W/D) cm	192 / 70 / 74.2
Gross / Net weight kg	149 / 135
Voltage V/Hz	220-240 / 50



Activated Charcoal Filter
Part no.: 988129100

* Noise output - according to EN 60704-3
 ** These models are only bound for storing wine. Max. DIN norm bordeaux bottle capacity (0.75 l) as per norm NF H 35-124 (H:300.5 mm, Ø 76.1mm)

Subject to modification. Current data see home.liebherr.com