



ENERGY STAR CERTIFIED

Lab Grade Refrigerators and Freezers

Liebherr - SUFsg 5001-70* : SUFsg 5001-70*

Specifications	
ENERGY STAR Unique ID:	2367456
ENERGY STAR Partner:	Liebherr-Hausgerate Ochsenhausen GmbH
Brand Name:	Liebherr
Model Name:	SUFsg 5001-70*
Model Number:	SUFsg 5001-70*
Product Type:	Ultra-Low Temperature Freezer
Total Volume (cubic feet):	17.7
Ultra-Low Temperature Freezer Energy Consumption (kWh/day/cu-ft):	0.46
Peak Temperature Variance (°C):	30.32
Product Form Factor:	Upright
Total Number of Doors:	1
Number of glass doors:	0
Number of solid doors:	1
Door Opening Orientation:	Swinging
Number of Drawers or Shelves:	3
Number of Internal Ports and/or Access Holes:	1
Depth (in):	39.0
Height:	77.5
Width:	35.0
Defrost Type:	Manual
Refrigerant with GWP:	R-170 (ethane) (GWP:6 Lower GWP), R-290 (GWP:3 Lower GWP)
Installed Accessories:	N/A
Date Available on Market:	1598932800000
Date Certified:	1602734400000
Markets:	United States, Canada
ENERGY STAR Certified:	Yes

Captured On:
06/10/2024