


LIEBHERR



Optimum humidity: Liebherr wine cabinets guarantee an optimum humidity and this ensured the corks are supple in place whilst still allowing the wine to breathe in its maturing process.




LED lighting: The tinted glass door provides not only a perfect UV protection, but also an optimum glance at wine. The integrated LED lighting ensures a pleasant illumination of the complete interior.

UV radiation. The LED lighting ensures that the interior is illuminated evenly and pleasantly so wines can be presented with illumination even over a longer period. WKb 1712

Barrique	
°Energy consumption (kWh/yr) I 24hr	°149 I 0.408
Bottle capacity (750ml Bordeaux)	63
Safety glass door	•
Control system Display	Buttons Digital
Digital temperature display	External
Tinted glazed door	Double
Temperature zones	1
Temperature range (°C) Min I Max °C	5 20
Suitable for ambient temperature Min I Max °C	10 43
Interior light	Vertical LED
Beech wood shelves	2
Base shelves	1
Recirculated air cooling Fresh air supply via activated charcoal filter	• 1
Child proof lock (control panel)	•
Handle design	Black stainless steel
Key lock Label holders	• •
Temp alarm I Door open warning signal	• •
Door hinge	Right, reversible
Gross weight (kg) I Net weight (kg)	50 45
Voltage (V) Rated Current (A)	220-240 1.0
Made In	Malaysia



Wooden shelves: The wooden shelves are ideal for storing Bordeaux bottles. They can be stored in opposite directions to make full use of the capacity.



 $^{\circ}$ = Measurements taken at +25°C ambient and +12°C interior temperature tested by Liebherr Germany (not to Australian standards).

Subject to modification. Current data see home.liebherr.com.au (240528)

Activated Charcoal Filter Part no.: YL9192811