

BLU  
Performance

# Legendary German Engineering

Quality, Design and Innovation



SWTNes 4265 & SKBes 4360



# LIEBHERR



**2.4" Touchscreen Display:** Located behind the appliance door, this high-resolution and high-contrast 2.4" touchscreen display can be read from any angle. Stored menu options are intuitive and easy-to-operate.



**BioFresh-Plus:** With the adjustable temperature range of 0°C and -2°C, BioFresh-Plus technology offers even more flexibility for your produce. The BioFresh-Plus drawer includes a versatile Fish & Seafood-Safe.





**PowerCooling:** A high-performance system that allows freshly stored items to be chilled rapidly while maintaining an even temperature throughout the interior. In this energy-efficient BluPerformance range, the fan switches off when the appliance door is opened, saving valuable energy.

*Note: Whilst every effort is made to ensure all specifications and information in this flyer are accurate, Liebherr reserves the right to make changes without notification. Where actual appliance measurements, specifications, and ratings are crucial, we recommend that the actual appliance installation instructions that accompany each appliance be checked. Liebherr will not be liable for any loss.*

Introducing Liebherr's award-winning Side-by-Side model from the BluPerformance range – the **SBSes 8486**. This eye-catching unit boasts an array of innovative solutions with six different climate zones, but its multi-temperature wine compartment is what makes it most unique! Enjoy added flexibility with BioFresh-Plus safe which offers an adjustable temperature range of 0°C and -2°C. Comes with a versatile Fish & Seafood-Safe.

<b>SBSes 8486</b> PremiumPlus	Bio Fresh-Plus	No Frost	Ice Maker	Wine	Soft System	Soft Telescopic	Smart Steel
----------------------------------	----------------	----------	-----------	------	-------------	-----------------	-------------

*Side-by-Sides comprised of two models	SWTNeS 4265	SKBes 4360
Energy consumption (kWh/yr)   Star rating	211   5.0	130   5.5
Compressors	1	1
Total Usable Litres	279	372
Refrigerator Usable Litres	-	236
Freezer Usable Litres	126	-
BioFresh Usable Litres	-	136
Wine Usable Litres	153 (48 bottles)	-
Control system and display   Child proof	2.4" Touch display   Yes	2.4" Touch display   Yes
PremiumPlus lighting	Yes	Yes
Sabbath mode	Yes	Yes
SoftSystem hinge - assists closure	Yes	Yes
Door open   Temperature : alarms	Yes   Yes	Yes   -
BioFreshPlus drawers on telescopic rails	-	1
BioFresh drawers on telescopic rails	-	2
SuperCool   VarioSafe	-   -	Yes   Yes
SuperFrost   VarioSpace	Yes   Yes	-   -
Integrated bottle shelf	-	Yes
Freezer drawers   Ice cube drawer	4   1	-   -
IceMaker   Ice cube capacity (kg/24h)	Yes   1.2	-   -
Stainless steel hose supplied (m)	1.5	-
Stainless steel with SmartSteel Doors   Walls	Yes   Yes	Yes   Yes
Default hinge	Left	Right
Voltage (V)   Rated Current (A)	220 - 240   1.4	220 - 240   1.2
Made in	Germany	Germany



\* These models are only bound for storing wine. Max. DIN norm bordeaux bottle capacity (0.75 l) as per norm NF H 35-124 (H:300.5 mm, Ø 76.1mm)

Important information: For Side-by-Side installation the freezer or the appliance with freezer compartment must be installed on the left as viewed from the front to prevent possible formation of condensation.