

CNef 4315

Fridge-freezer with BioCool and NoFrost
Comfort



Freestanding

CNef 4315

Fridge-freezer with BioCool and NoFrost
Comfort

← 60 →



Energy consump. in 365 days/24 h	220 / 0.602 kWh ¹
Volume	Freezer compartment: 120 l
Climate rating	SN-T (+10 °C to +43 °C)
Frequency / voltage	50/60 Hz / 220-240 V~
Connection Rating	1.4 A 190 W
Door / Side walls	SmartSteel / silver
Exterior dimensions in cm (h / w / d)	185 / 60 / 66.5 ¹
Controls	
Type of Control	2.4" TFT-monochrome display behind-the-door, touch electronic
Interface	SmartDeviceBox: optional accessory
Door alarm	visual and acoustic
Alarm in case of malfunction	visual and acoustic
Fridge compartment	
Forced air system	PowerCooling with FreshAir activated charcoal filter
Door interior	Comfort GlassLine storage shelf with bottle holder, Egg tray
Bottle shelf	Bottle shelf
Interior light	LED ceiling lighting
Fruit and vegetable bins	1 BioCool-Box on roller tracks
Freezer compartment	
Star coding	****
Number of drawers	3 drawers, VarioSpace
Temperature rise time	24 h
Ice cube tray	yes
Features	
Door hinges	Reversible door hinging
EAN-Nr.	4016803050612

¹ To achieve the indicated energy consumption level, the wall spacers provided with the appliance must be used. These increase the appliance depth by approximately 1.5 cm. If the wall spacers are not used, the appliance will remain fully functional, but will have a slightly higher energy consumption level.

Features



BluPerformance

Offering even more interior capacity for food storage, economy and energy efficiency, quiet operation, ease-of-use and exceptional design: Liebherr's BluPerformance appliances set both a new direction in refrigeration and freezing, and new benchmarks.



BioCool

The BioCool-Box is a feature of the refrigerator compartment and the humidity level within it can be regulated to keep fruit and vegetables fresher for longer. It extends smoothly on either a roller system or on premium-quality telescopic rails, both of which guarantee ease-of-use and allow convenient access for a perfect view of stored products. Humidity is regulated via slider control.



NoFrost

The NoFrost system from Liebherr provides rapid, fan-assisted, frost-free freezing for the safe, longterm storage of food. The NoFrost technology creates a much larger storage capacity and keeps the freezer constantly frost-free.




DuoCooling

DuoCooling uses two completely separate refrigeration circuits to ensure no air is exchanged between the fridge and freezer compartments. The food neither dries out nor is odour transferred. This means throwing less away, shopping less often, but saving more and greater enjoyment.



VarioSpace


Freezer drawers and toughened glass shelves can be removed. VarioSpace – the practical storage solution for larger items.



SuperFrost

SuperFrost generates greater cold reserves so that freshly stored food can be frozen rapidly, in a vitamin-preserving manner. Reversion from the freezing temperature of -32 °C to the normal freezer temperature is time or quantity-controlled to help save energy.

Technical drawing



Products and their details are staged for advertising purposes. Whilst every effort is made to ensure all specifications and information in this flyer are accurate, Liebherr reserves the right to make changes without notification. Where actual appliance measurements, specifications, and ratings are crucial, we recommend that the actual appliance installation instructions that accompany each appliance be checked. Liebherr will not be liable for any damage. (01.12.2022)

