




ELEVATION LEFT UWTes 1672
1 : 10



ELEVATION FRONT UWTes 1672
1 : 10



FLOOR PLAN UWTes 1672
1 : 5



3D View UWTes 1672

UWTes 1672 Vinidor	
Total capacity litre	94
Capacity of wine compartment litre	94
Max. No. of Bordeaux btls 0.75	34 *
Shipping weight kg	57
Climate rating	T
Door hinges	right, reversible
Outside dimensions in cm (H / W / D)	82-87 / 60 / at least 58
Energy supply Voltage V/Hz	220 – 240 / 50
Electronic control	MagicEye in the interior, key electronic Digital temperature display for both wine zones

* These models are only bound for storing wine. Maximum DIN norm bordeaux bottles capacity (0.75 l) as per norm NF H 35-124 (H: 300.5 mm, Ø 76.1 mm)

UWTes 1672_136

ScaleAs indicated