



ENERGY STAR CERTIFIED

Lab Grade Refrigerators and Freezers

Liebherr - SRFvh 5501 : SRFvh 5501

Specifications

ENERGY STAR Unique ID:	2668047
ENERGY STAR Partner:	Liebherr-Hausgerate Ochsenhausen GmbH
Brand Name:	Liebherr
Model Name:	SRFvh 5501
Model Number:	SRFvh 5501
Product Type:	High Performance Refrigerator
Total Volume (cubic feet):	18.37
Lab Grade Refrigerator or Freezer Energy Consumption (kWh/day):	0.77
Peak Temperature Variance (°C):	5.36
Product Form Factor:	Upright
Total Number of Doors:	1
Number of glass doors:	0
Number of solid doors:	1
Door Opening Orientation:	Swinging
Number of Drawers or Shelves:	5
Number of Internal Ports and/or Access Holes:	1
Depth (in):	30.71
Height:	70.55
Width:	29.41
Defrost Type:	Automatic
Refrigerant with GWP:	R-600a (GWP:3)
Installed Accessories:	n/a
Date Available on Market:	1706763600000
Date Certified:	1695182400000
Markets:	United States, Canada
ENERGY STAR Certified:	Yes

Captured On:
06/10/2024