

American style fridge-freezers

Freestanding

Model	+ 3 years	+ 8 years
SBSbs8673	£ 225	£ 400
SBSbs8683	£ 225	£ 400
SBSef7242	£ 190	£ 360
SBSes7165	£ 225	£ 400
SBSes7353	£ 190	£ 360
SBSes8473	£ 225	£ 400
SBSes8483	£ 225	£ 400
SBSes8486	£ 225	£ 400
SBSes8496	£ 225	£ 400
SBSes8663	£ 225	£ 400
SBSes8773	£ 225	£ 400

Built-In

Model	+ 3 years	+ 8 years
IXCC5155	£ 225	£ 400
IXCC5165	£ 225	£ 400
SBS70I2	£ 225	£ 400
SBS70I4	£ 225	£ 400

Fridge-freezers

Freestanding

Model	+ 3 years	+ 8 years
CBN4815	£ 170	£ 330
CBNbe6256	£ 225	£ 400
CBNbs4815	£ 170	£ 330
CBNef4815	£ 170	£ 330
CBNef5715	£ 170	£ 330
CBNef5735	£ 170	£ 330
CBNes6256	£ 225	£ 400
CBNesf5133	£ 145	£ 275
CBNP4858	£ 170	£ 330
CBNPbs4858	£ 170	£ 330
CBNPes3967	£ 170	£ 330
CBNPes4858	£ 170	£ 330
CBNPes4878	£ 170	£ 330
CBNPes5156	£ 170	£ 330
CBNPes5167	£ 170	£ 330
CBNPes5758	£ 170	£ 330
CBNPgb4855	£ 170	£ 330
CBNPgw4855	£ 170	£ 330

Built-In

Model	+ 3 years	+ 8 years
ECBN9471	£ 225	£ 400
ECBN9671	£ 225	£ 400
ECBN9673	£ 225	£ 400
ECBN5066	£ 225	£ 400
ECBN6156	£ 225	£ 400
ECBN6256	£ 225	£ 400
ICBb5152	£ 145	£ 275
ICBd5122	£ 125	£ 240
ICBdi5122	£ 125	£ 240
ICBN3324	£ 145	£ 275
ICBN3376	£ 170	£ 330
ICBN3386	£ 170	£ 330
ICBNd5153	£ 145	£ 275
ICBNd5163	£ 145	£ 275
ICBNdi5183	£ 170	£ 330
ICBNe5123	£ 125	£ 240
ICBP3266	£ 170	£ 330
ICd5123	£ 125	£ 240
ICe5103	£ 60	£ 120
ICN3314	£ 125	£ 240
ICN3386	£ 170	£ 330
ICNd5123	£ 125	£ 240

Fridge-freezers

Freestanding

Model	+ 3 years	+ 8 years
CN3033	£ 125	£ 240
CN3115	£ 125	£ 240
CN3515	£ 125	£ 240
CN3915	£ 125	£ 240
CN4015	£ 145	£ 275
CN4213	£ 60	£ 120
CN4313	£ 60	£ 120
CN4315	£ 145	£ 275
CN4713	£ 60	£ 120
CN4813	£ 60	£ 120
CN4815	£ 145	£ 275
CN5715	£ 145	£ 275
CNbe4015	£ 145	£ 275
CNbs3915	£ 145	£ 275
CNbs4015	£ 145	£ 275
CNbs4315	£ 170	£ 330
CNef3115	£ 125	£ 240
CNef3515	£ 145	£ 275
CNef3535	£ 145	£ 275
CNef3915	£ 145	£ 275
CNef4015	£ 145	£ 275
CNef4315	£ 145	£ 275
CNef4815	£ 145	£ 275
CNef5715	£ 170	£ 330
CNef5725	£ 145	£ 275
CNef5735	£ 145	£ 275
CNef5745	£ 145	£ 275
CNel4213	£ 60	£ 120
CNel4313	£ 60	£ 120
CNel4713	£ 60	£ 120
CNel4813	£ 60	£ 120
CNfb4313	£ 60	£ 120
CNkw4313	£ 60	£ 120
CNno4313	£ 60	£ 120
CNP4313	£ 60	£ 120
CNP4758	£ 170	£ 330
CNP4813	£ 60	£ 120
CNPel4313	£ 60	£ 120
CNPel4813	£ 60	£ 120
CNPes4758	£ 170	£ 330

Built-In

Model	+ 3 years	+ 8 years
ICNdi5173	£ 170	£ 330
ICNe5133	£ 125	£ 240
ICNf5103	£ 60	£ 120
ICNS3324	£ 125	£ 240
ICNSf5103	£ 60	£ 120
ICP3324	£ 125	£ 240
ICSe5103	£ 60	£ 120
ICSe5122	£ 125	£ 240
ICU3324	£ 125	£ 240
ICUN3324	£ 125	£ 240
ICUNS3324	£ 125	£ 240
ICUS2924	£ 125	£ 240
ICUS3224	£ 125	£ 240
ICUS3324	£ 125	£ 240
SICN3386	£ 145	£ 275
SICNd5153	£ 145	£ 275

Fridge-freezers

Freestanding

Model	+ 3 years	+ 8 years
CNPes4858	£ 170	£ 330
CNPes4868	£ 170	£ 330
CNsl3033	£ 125	£ 240
CP4313	£ 60	£ 120
CP4813	£ 60	£ 120
CPel4313	£ 60	£ 120
CPel4813	£ 60	£ 120
CU2311	£ 60	£ 120
CU2331	£ 60	£ 120
CU2811	£ 60	£ 120
CU2831	£ 60	£ 120
CU3311	£ 60	£ 120
CU3331	£ 60	£ 120
CUag3311	£ 60	£ 120
CUef2811	£ 60	£ 120
CUef3311	£ 60	£ 120
CUel2331	£ 60	£ 120
CUel2831	£ 60	£ 120
CUel3331	£ 60	£ 120
CUfb2831	£ 60	£ 120
CUfr3311	£ 60	£ 120
CUkw2831	£ 60	£ 120
CUN3033	£ 125	£ 240
CUno2831	£ 60	£ 120
CUsl2311	£ 60	£ 120
CUsl2811	£ 60	£ 120
CUwb3311	£ 60	£ 120

Top mount freezer

Model	+ 3 years	+ 8 years
CT2131	£ 60	£ 120
CT2531	£ 60	£ 120
CT2931	£ 60	£ 120
CTel2131	£ 60	£ 120
CTel2531	£ 60	£ 120
CTel2931	£ 60	£ 120
CTN3663	£ 145	£ 275
CTN5215	£ 145	£ 275
CTNef5215	£ 145	£ 275
CTNesf3663	£ 145	£ 275
CTP2121	£ 60	£ 120
CTP2521	£ 60	£ 120
CTP2921	£ 60	£ 120
CTP3316	£ 125	£ 240
CTPag2121	£ 60	£ 120
CTPesf3316	£ 145	£ 275
CTPfr2121	£ 60	£ 120
CTPsl2121	£ 60	£ 120
CTPsl2521	£ 60	£ 120
CTPwb2121	£ 60	£ 120

Tall fridges

Freestanding

Model	+ 3 years	+ 8 years
K2330	£ 70	£ 140
K2340	£ 70	£ 140
K2630	£ 80	£ 165
K2814	£ 70	£ 140
K2834	£ 70	£ 140
K3130	£ 80	£ 165
K3710	£ 100	£ 190
K3730	£ 100	£ 190
K4310	£ 100	£ 190
K4330	£ 100	£ 190
KB3750	£ 115	£ 225
KB4310	£ 115	£ 225
KB4330	£ 115	£ 225
KB4350	£ 115	£ 225
KBbs4350	£ 115	£ 225
KBef4310	£ 115	£ 225
KBef4330	£ 115	£ 225
KBes3750	£ 115	£ 225
KBes4350	£ 115	£ 225
KBPes4354	£ 115	£ 225
KBPgb4354	£ 115	£ 225
KBPgw4354	£ 115	£ 225
Kef3710	£ 115	£ 225
Kef3730	£ 115	£ 225
Kef4310	£ 115	£ 225
Kef4330	£ 115	£ 225
Ksl2630	£ 80	£ 165
Ksl2814	£ 70	£ 140
Ksl2834	£ 70	£ 140
Ksl3130	£ 80	£ 165

Built In

Model	+ 3 years	+ 8 years
EKB9271	£ 225	£ 400
EKB9471	£ 225	£ 400
EKB9671	£ 225	£ 400
IK1960	£ 100	£ 190
IK2320	£ 80	£ 165
IK2324	£ 80	£ 165
IK2360	£ 100	£ 190
IK2760	£ 100	£ 190
IK2764	£ 100	£ 190
IK3520	£ 80	£ 165
IK3524	£ 80	£ 165
IKB2360	£ 115	£ 225
IKB3520	£ 100	£ 190
IKB3560	£ 115	£ 225
IKBP2760	£ 115	£ 225
IKBP2964	£ 115	£ 225
IKBP3560	£ 115	£ 225
IKBP3564	£ 115	£ 225
IKF3510	£ 100	£ 190
IRBb4170	£ 115	£ 225
IRBd4020	£ 80	£ 165
IRBd4050	£ 100	£ 190
IRBd4150	£ 100	£ 190
IRBd4550	£ 100	£ 190
IRBd4570	£ 115	£ 225
IRBd5151	£ 100	£ 190
IRBdi5180	£ 115	£ 225
IRBe5120	£ 80	£ 165
IRBSe5120	£ 80	£ 165
IRd3950	£ 100	£ 190
IRe4100	£ 70	£ 140
IRe4101	£ 70	£ 140
IRe4520	£ 80	£ 165
IRe4521	£ 80	£ 165
IRf3900	£ 70	£ 140
IRf5101	£ 70	£ 140
IRSF3900	£ 70	£ 140
IRSF3901	£ 70	£ 140

Under counter fridges

Freestanding

Model	+ 3 years	+ 8 years
OKes1750	£ 145	£ 275
T1400	£ 45	£ 80
T1404	£ 45	£ 80
T1410	£ 55	£ 95
T1414	£ 55	£ 95
T1504	£ 45	£ 80
T1700	£ 45	£ 80
T1710	£ 55	£ 95
T1714	£ 55	£ 95
T1810	£ 55	£ 95
Tb1400	£ 45	£ 80
TP1410	£ 55	£ 95
TP1414	£ 55	£ 95
TP1424	£ 55	£ 95
TP1434	£ 55	£ 95
TP1720	£ 55	£ 95
TP1724	£ 55	£ 95
TP1760	£ 65	£ 110
TP1764	£ 65	£ 110
TPesf1710	£ 65	£ 110
TPesf1714	£ 65	£ 110
Tsl1414	£ 55	£ 95
TX1021	£ 55	£ 95

Built In

Model	+ 3 years	+ 8 years
BP2850	£ 115	£ 225
EK1620	£ 80	£ 165
IB1650	£ 100	£ 190
IK1620	£ 80	£ 165
IK1624	£ 80	£ 165
IKP1660	£ 100	£ 190
IKS1620	£ 80	£ 165
IKS1624	£ 80	£ 165
SIBP1650	£ 100	£ 190
SUIB1550	£ 115	£ 225
UIK1424	£ 80	£ 165
UIK1510	£ 80	£ 165
UIK1514	£ 80	£ 165
UIK1550	£ 100	£ 190
UIK1610	£ 80	£ 165
UIK1620	£ 80	£ 165
UIKo1550	£ 115	£ 225
UIKo1560	£ 115	£ 225
UIKP1550	£ 100	£ 190
UIKP1554	£ 100	£ 190
UK1524	£ 80	£ 165
UK1720	£ 80	£ 165
UKes1752	£ 145	£ 275

Tall freezers

Freestanding

Model	+ 3 years	+ 8 years
G3513	£ 70	£ 140
G4013	£ 70	£ 140
GN2835	£ 100	£ 190
GN3235	£ 80	£ 165
GN3613	£ 80	£ 165
GN3735	£ 100	£ 190
GN4113	£ 80	£ 165
GN4335	£ 100	£ 190
GN4375	£ 100	£ 190
GN5215	£ 100	£ 190
GN5275	£ 100	£ 190
GNef4335	£ 100	£ 190
GNP1913	£ 80	£ 165
GNP2313	£ 80	£ 165
GNP2713	£ 80	£ 165
GNP3013	£ 80	£ 165
GNP3755	£ 115	£ 225
GNP4355	£ 115	£ 225
GNPes4355	£ 115	£ 225
GP2033	£ 70	£ 140
GP2433	£ 70	£ 140
GP2733	£ 70	£ 140
GP3013	£ 80	£ 165
GP3513	£ 80	£ 165
GP4013	£ 80	£ 165

Built In

Model	+ 3 years	+ 8 years
EGN9171	£ 225	£ 400
EGN9271	£ 225	£ 400
EGN9471	£ 225	£ 400
EGN9671	£ 225	£ 400
SIGN2756	£ 115	£ 225
SIGN3524	£ 100	£ 190
SIGN3556	£ 115	£ 225
SIGN3576	£ 115	£ 225
SIFNe5188	£ 115	£ 225

Under counter freezers

Freestanding

Model	+ 3 years	+ 8 years
G1213	£ 45	£ 80
G1223	£ 45	£ 80
GN1066	£ 65	£ 110
GNP1066	£ 65	£ 110
GNsl2323	£ 80	£ 165
GP1213	£ 45	£ 80
GP1376	£ 55	£ 95
GP1476	£ 55	£ 95
GP1486	£ 65	£ 110
GPesf1476	£ 65	£ 110
Gsl1223	£ 45	£ 80
GX823	£ 55	£ 95

Built In

Model	+ 3 years	+ 8 years
IG1024	£ 80	£ 165
IG1624	£ 80	£ 165
IGN1064	£ 100	£ 190
IGN1664	£ 100	£ 190
IGS1624	£ 80	£ 165
SUIG1514	£ 80	£ 165
SUIGN1554	£ 100	£ 190
UIG1313	£ 80	£ 165
UIG1323	£ 80	£ 165

Chest freezers

Model	+ 3 years	+ 8 years
GT1432	£ 80	£ 165
GT2132	£ 80	£ 165
GT2632	£ 80	£ 165
GT3632	£ 80	£ 165
GT6122	£ 80	£ 165

Tall wine cabinets

Freestanding

Model	+ 3 years	+ 8 years
WKb3212	£ 145	£ 275
WKb4212	£ 145	£ 275
WKes4552	£ 170	£ 330
WKr3211	£ 145	£ 275
WKr4211	£ 145	£ 275
Wkt4551	£ 170	£ 330
Wkt4552	£ 170	£ 330
Wkt5551	£ 170	£ 330
Wkt5552	£ 170	£ 330
Wkt6451	£ 170	£ 330
WTb4212	£ 170	£ 330
WTes5872	£ 170	£ 330
WTes5972	£ 170	£ 330
WTpes5972	£ 170	£ 330
WTr4211	£ 145	£ 275

Built In

Model	+ 3 years	+ 8 years
EWT9175	£ 225	£ 400
EWT9275	£ 225	£ 400
EWTdf1653	£ 170	£ 330
EWTdf2353	£ 170	£ 330
EWTdf3553	£ 170	£ 330
EWTgb1683	£ 170	£ 330
EWTgb2383	£ 170	£ 330
EWTgb3583	£ 170	£ 330
EWTgw1683	£ 170	£ 330
EWTgw2383	£ 170	£ 330
EWTgw3583	£ 170	£ 330

Under counter wine cabinets

Freestanding

Model	+ 3 years	+ 8 years
WKb1812	£ 145	£ 275
WTes1672	£ 170	£ 330

Built In

Model	+ 3 years	+ 8 years
UWKes1752	£ 145	£ 275
UWT1682	£ 170	£ 330
UWTes1672	£ 170	£ 330
UWTgb1682	£ 170	£ 330
WKEes553	£ 170	£ 330
WKEgb582	£ 170	£ 330
WKEgw582	£ 170	£ 330
WKes653	£ 170	£ 330

Humidor

Model	+ 3 years	+ 8 years
ZKes453	£ 170	£ 330