

Liebherr - freshness power with BioFresh App

- The iOS and Android-compatible BioFresh app for Smart phones and tablets provides first class service and comprehensive advice
- BioFresh by Liebherr is the perfect technology for ultra-long enjoyment of fresh food

August 2013 – Appliances featuring BioFresh technology are becoming increasingly popular as they extend the storage life of food and a noticeable quality of freshness is experienced every day. BioFresh technology strongly appeals to customers by providing both, day-to-day convenience and impressively long-lasting freshness.

The iOS and Android-compatible BioFresh app for Smart phones and tablets provides first class service and comprehensive advice. Liebherr's new, free 'BioFresh HD' app for tablets is simple and intuitive to use and provides information about shelf life and storage, as well as about the valuable vitamins and minerals contained within foods and their contribution to everyday health and fitness.

Download the Liebherr BioFresh app and find out lots of valuable information about your food, its vitamins and shelf life. Search for your favourite food or products that contain a specific vitamin – in alphabetic order or by category.

Captions

Image: Liebherr_BioFresh App_Icon

Liebherr's new, free 'BioFresh HD' app for tablets is simple and intuitive to use and provides information about shelf life and storage, as well as about the valuable vitamins and minerals contained within foods and their contribution to everyday health and fitness.

Contact person

Claudia Meindl

Marketing & Communication

Phone: +43 2262 602-454

E-mail: claudia.meindl@liebherr.com

Published by

Liebherr-Hausgeräte GmbH

Ochsenhausen / Deutschland

www.liebherr.com